

GrpE E.Coli

Description: GrpE Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 197 amino acids and having a molecular mass of 21.8 kDa. The GrpE is purified by proprietary chromatographic techniques.

Catalog #:HYP5-021

For research use only.

Synonyms: HSP24, HSPB25.3, HSP-70 Cofactor, Protein grpE, Heat shock protein B25.3, grpE, b2614, JW2594.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MSSKEQKTPE GQAPEEIIMD QHEEIEAVEP EASAEQVDPR
DEKIANLEAQ LAEAQTRERD GILRVKAEME NLRRTRELDI EKAHKFALEK FINELLPVID
SLDRALEVAD KANPDMSAMV EGIELTLKSM LDVVRKFGVE VIAETNVPLD PNVHQAIAMV
ESDDVAPGNV LGIMQKGYTL NGRTIRAAMV TVAKAKA.

Purity: Greater than 90.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

The HSP24 protein solution contains 20mM Tris-HCl pH-8 and 0.1M NaCl.

Stability:

GrpE although stable 4°C for 4 weeks, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

GrpE, co-chaperone of E.coli, participates actively in the response to hyperosmotic and heat shock by preventing the aggregation of stress-denatured proteins in association with DnaK. This protein is the nucleotide exchange factor for DnaK and may function as a thermosensor. Several rounds of ATP-dependent interactions between DnaJ, DnaK and GrpE are required for fully efficient folding.

To place an order, please [Click HERE](#).